


Antigen binding-avidity antibody titer assay

AMED Support Project JP19he2302005

[Advantages of the Novel Inspection]

- Data predict oral food challenge (OFC) -positivity or -negativity in **individual patients**.
- Detection of the low binding-avidity IgE antibodies and OFC-negativity, **providing future safe dietary intake doses in dietary life**.
- Monitor individual **oral immunotherapy efficacy** over time.


1. What is antigen binding-avidity antibody titer (ABAT) ?



ABAT (antigen binding-avidity antibody titer)

Evaluation by **Quantity x Quality** instead of conventional antibody quantity measurement.



2. Measurement of binding-avidity of antibodies by competitive binding inhibition between immobilized antigen on the DCP microarray and added serially diluted soluble antigen.



Antigen binding-avidity = 1/IC₅₀ (nM)

3. Identification of OFC negativity in allergic patients

Nutrients, 2023,15(12):2770



In the figure on the right, which shows sIgE/IC₅₀ as a reaction sensitivity and severity biomarker and the ratio of sIgE/sIgG4 as a reaction onset biomarker on the Y-axis and X-axis, respectively. Each cut-off values based on ROC analysis of OFC negativity/positivity determine individual OFC outcomes. **Negative probability in the D-area provides future safe dietary intake doses in daily life.**

Although the current ImmunoCAP OVM sIgE assay detects statistically significant differences between OFC positive- and negative-groups, **it is difficult to determine the OFC-positivity or -negativity of individual patients.**

Allergen-specific antibody assay (IgG1, IgG4, IgA, IgE)
Antigen binding-avidity antibody titer assay

for Reserch use

	Allergen				Sample Volume	Storage
Allergen-specific antibody assay	Egg White	House dust	Ovomucoid	Der p 1	Serum	-25°C
	Egg York	Dermatophagoides farina	Ovalbumin	Der p 2		
Specific IgE	Cow Milk	Dermatophagoides pteronyssinus	Casein	Alt a 1	100 μL	-25°C
	Buckwheat	Dog dander	α-Casein	Cl a h 8		
Specific IgG1	Peanut	Cat dander	β-Casein	Asp r 1	100 μL	-25°C
	Shrimp	Japanese Cedar	β-Lactoglobulin			
Specific IgG4	Crab	Hinoki Cypress	Gluten		100 μL	-25°C
	Salmon	Birch White	Gliadin			
Specific IgA	Tuna	Alder European	Ara h 1		100 μL	-25°C
		Anthoxanthum odoratum	Ara h 2			
Antigen binding-avidity antibody titer assay		Orchard grass	Cry j 1		100 μL	-25°C
		Ragweed	Cry j 2			
		Timothy	Der f 1		100 μL	-25°C
			Der f 2			



Contract manufacturing and inspection of DCP chips

For inquiries or requests for inspections, please contact us at the contact information below.

Contact us :
Applied Medical Enzyme Research Institute Corporation

U R L : <http://www.americ.co.jp/>